

## Pre-mRNA-splicing regulator WTAP

Catalog No: #AP82736



Package Size: #AP82736-1 50ug #AP82736-2 100ug #AP82736-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Pre-mRNA-splicing regulator WTAP  |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Hu  |
| Other Names        | Female-lethal(2)D homolog,hFL(2)D,WT1-associated protein,Wilms tumor 1-associating protein,KIAA0105 |
| Accession No.      | Uniprot ID:Q15007   |
| Uniprot            | Q15007  |
| GeneID             | 9589;   |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)  |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:WTAP

Note: This product is for in vitro research use only