

Filamin-A

Catalog No: #AP82737



Package Size: #AP82737-1 50ug #AP82737-2 100ug #AP82737-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Filamin-A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FLN-A,Actin-binding protein 280,ABP-280,Alpha-filamin,Endothelial actin-binding protein,Filamin-1,Non-muscle filamin,Fln,Fln1
Accession No.	Uniprot ID:Q8BTM8
Uniprot	Q8BTM8
GeneID	192176;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Flna

Note: This product is for in vitro research use only