

Ephrin-A3

Catalog No: #AP82740



Package Size: #AP82740-1 50ug #AP82740-2 100ug #AP82740-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Ephrin-A3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | EHK1 ligand,EHK1-L,EPH-related receptor tyrosine kinase ligand 3,LERK-3,Epl3,Eplg3,Lerk3 |
| Accession No. | Uniprot ID:O08545 |
| Uniprot | O08545 |
| GeneID | 13638; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Efna3

Note: This product is for in vitro research use only