

## DNA polymerase alpha subunit B

Catalog No: #AP82741



Package Size: #AP82741-1 50ug #AP82741-2 100ug #AP82741-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                        |
|--------------------|--------------------------------------------------------|
| Product Name       | DNA polymerase alpha subunit B                         |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli                                                 |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Hu                                                     |
| Other Names        | DNA polymerase alpha 70 kDa subunit                    |
| Accession No.      | Uniprot ID:Q14181                                      |
| Uniprot            | Q14181                                                 |
| GeneID             | 23649;                                                 |
| Tag Info           | His                                                    |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: POLA2

Note: This product is for in vitro research use only