

EF-hand domain-containing protein D2

Catalog No: #AP82745



Package Size: #AP82745-1 50ug #AP82745-2 100ug #AP82745-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	EF-hand domain-containing protein D2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Swiprosin-1,SWS1
Accession No.	Uniprot ID:Q96C19
Uniprot	Q96C19
GeneID	79180;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EFHD2

Note: This product is for in vitro research use only