

## Fatty acid-binding protein, brain

Catalog No: #AP82757



Package Size: #AP82757-1 50ug #AP82757-2 100ug #AP82757-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fatty acid-binding protein, brain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Brain-type fatty acid-binding protein,B-FABP,Fatty acid-binding protein 7,Fatty acid-binding protein, retina,R-FABP
Accession No.	Uniprot ID:Q05423
Uniprot	Q05423
GeneID	396246;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FABP7

Note: This product is for in vitro research use only