

Glyceraldehyde-3-phosphate dehydrogenase

Catalog No: #AP82758



Package Size: #AP82758-1 50ug #AP82758-2 100ug #AP82758-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glyceraldehyde-3-phosphate dehydrogenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GAPDH,38 kDa BFA-dependent ADP-ribosylation substrate,BARS-38,Peptidyl-cysteine S-nitrosylase GAPDH,Gapd
Accession No.	Uniprot ID:P04797
Uniprot	P04797
GeneID	108351137;24383;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gapdh

Note: This product is for in vitro research use only