

Eosinophil cationic protein

Catalog No: #AP82762



Package Size: #AP82762-1 50ug #AP82762-2 100ug #AP82762-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Eosinophil cationic protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ECP,Eosinophil secondary granule ribonuclease 11,EAR-11,Ribonuclease 3,RNase 3,Ear11,Rns3
Accession No.	Uniprot ID:P70709
Uniprot	P70709
GeneID	100361866;192264;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rnase3

Note: This product is for in vitro research use only