

Pikachurin

Catalog No: #AP82763



Package Size: #AP82763-1 50ug #AP82763-2 100ug #AP82763-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pikachurin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	EGF-like, fibronectin type-III and laminin G-like domain-containing protein
Accession No.	Uniprot ID:A3KN33
Uniprot	A3KN33
GeneID	534427;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: EGFLAM

Note: This product is for in vitro research use only