

# EGF, latrophilin and seven transmembrane domain-containing protein 1

Catalog No: #AP82764

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP82764-1 50ug #AP82764-2 100ug #AP82764-3 1mg

## Description

Product Name	EGF, latrophilin and seven transmembrane domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EGF-TM7-latrophilin-related protein,ETL protein,ETL,UNQ202/PRO228
Accession No.	Uniprot ID:Q9HBW9
Uniprot	Q9HBW9
GeneID	64123;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ELTD1

Note: This product is for in vitro research use only