

## Ephrin type-B receptor 4

Catalog No: #AP82766



Package Size: #AP82766-1 50ug #AP82766-2 100ug #AP82766-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ephrin type-B receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Hepatoma transmembrane kinase,Tyrosine-protein kinase TYRO11,HTK,MYK1,TYRO11
Accession No.	Uniprot ID:P54760
Uniprot	P54760
GeneID	2050;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:EPHB4

Note: This product is for in vitro research use only