

Excitatory amino acid transporter 1

Catalog No: #AP82767



Package Size: #AP82767-1 50ug #AP82767-2 100ug #AP82767-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Excitatory amino acid transporter 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Sodium-dependent glutamate/aspartate transporter 1, GLAST-1, Solute carrier family 1 member 3, EAAT1, GLAST, GLAST1
Accession No.	Uniprot ID:P46411
Uniprot	P46411
GeneID	282354;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SLC1A3

Note: This product is for in vitro research use only