

Extracellular matrix protein 1

Catalog No: #AP82769

Package Size: #AP82769-1 50ug #AP82769-2 100ug #AP82769-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Extracellular matrix protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Secretory component p85
Accession No.	Uniprot ID:Q62894
Uniprot	Q62894
GeneID	116662;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ecm1

Note: This product is for in vitro research use only