

## Ephrin-A2

Catalog No: #AP82770



Package Size: #AP82770-1 50ug #AP82770-2 100ug #AP82770-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ephrin-A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CEK7-ligand,CEK7-L,ELF-1,EPH-related receptor tyrosine kinase ligand 6,LERK-6,Elf1,Epl6,Eplg6,Lerk6
Accession No.	Uniprot ID:P52801
Uniprot	P52801
GeneID	13637;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Efna2

Note: This product is for in vitro research use only