

Endothelin-2

Catalog No: #AP82773



Package Size: #AP82773-1 50ug #AP82773-2 100ug #AP82773-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Endothelin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ET-2,Preproendothelin-2,PPET2
Accession No.	Uniprot ID:Q867A9
Uniprot	Q867A9
GeneID	319094;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EDN2

Note: This product is for in vitro research use only