

Ephrin-A2

Catalog No: #AP82774



Package Size: #AP82774-1 50ug #AP82774-2 100ug #AP82774-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ephrin-A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EPH-related receptor tyrosine kinase ligand 6, LERK-6, HEK7 ligand, HEK7-L, EPLG6, LERK6
Accession No.	Uniprot ID:O43921
Uniprot	O43921
GeneID	1943;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: EFNA2

Note: This product is for in vitro research use only