

Adropin

Catalog No: #AP82783



Package Size: #AP82783-1 50ug #AP82783-2 100ug #AP82783-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Adropin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Energy homeostasis-associated protein
Accession No.	Uniprot ID:Q8K1D8
Uniprot	Q8K1D8
GeneID	69638;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Enho

Note: This product is for in vitro research use only