

EGF-like module-containing mucin-like hormone receptor-like 1

Catalog No: #AP82786

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP82786-1 50ug #AP82786-2 100ug #AP82786-3 1mg

Description

Product Name	EGF-like module-containing mucin-like hormone receptor-like 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cell surface glycoprotein F4/80,EGF-like module receptor 1,EMR1 hormone receptor,Gpf480
Accession No.	Uniprot ID:Q61549
Uniprot	Q61549
GeneID	13733;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Emr1

Note: This product is for in vitro research use only