

ASK1 (Phospho-Thr838) Antibody

Catalog No: #12463

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ASK1 (Phospho-Thr838) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms
Specificity	ASK1 (Phospho-Thr838) Antibody detects endogenous levels of ASK1 only when phosphorylated at Thr838
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ASK1 (Phospho-Thr838)
Target Name	ASK1
Modification	Phospho
Other Names	MAP3K5, ASK-1, ASK1, MEKK5, MEK kinase 5, MEKK 5, MAPK/ERK kinase kinase 5, MAPKKK5, MAP/ERK kinase kinase 5
Accession No.	Swiss-Prot#: Q99683NCBI Gene ID: 4217
Uniprot	Q99683
GeneID	4217;
Target Species	human
Calculated MW	155kd
Concentration	1.0mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of ASK1 (pT838) expression in K562 (A) whole cell lysates.

Note: This product is for in vitro research use only