

Disks large homolog 4

Catalog No: #AP82791



Package Size: #AP82791-1 50ug #AP82791-2 100ug #AP82791-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Disks large homolog 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Postsynaptic density protein 95, PSD-95, Synapse-associated protein 90, SAP-90, Dlgh4, Psd95
Accession No.	Uniprot ID:P31016
Uniprot	P31016
GeneID	29495;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Dlgh4

Note: This product is for in vitro research use only