

Ezrin

Catalog No: #AP82794



Package Size: #AP82794-1 50ug #AP82794-2 100ug #AP82794-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ezrin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Cytovillin,Villin-2,p81,VIL2
Accession No.	Uniprot ID:P31976
Uniprot	P31976
GeneID	281574;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EZR

Note: This product is for in vitro research use only