

Deoxyuridine 5'-triphosphate nucleotidohydrolase

Catalog No: #AP82798



Package Size: #AP82798-1 50ug #AP82798-2 100ug #AP82798-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Deoxyuridine 5'-triphosphate nucleotidohydrolase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	dUTPase,PPAR-interacting protein 4,PIP4,dUTP pyrophosphatase,Dutp
Accession No.	Uniprot ID:P70583
Uniprot	P70583
GeneID	497778;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dut

Note: This product is for in vitro research use only