

Eukaryotic translation initiation factor 2-alpha kinase 3

Catalog No: #AP82799



Package Size: #AP82799-1 50ug #AP82799-2 100ug #AP82799-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Eukaryotic translation initiation factor 2-alpha kinase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PRKR-like endoplasmic reticulum kinase,Pancreatic eIF2-alpha kinase,Pek,Perk
Accession No.	Uniprot ID:Q9Z1Z1
Uniprot	Q9Z1Z1
GeneID	29702;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Eif2ak3

Note: This product is for in vitro research use only