

Dual specificity protein phosphatase 4

Catalog No: #AP82803



Package Size: #AP82803-1 50ug #AP82803-2 100ug #AP82803-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dual specificity protein phosphatase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q8BFV3
Uniprot	Q8BFV3
GeneID	319520;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dusp4

Note: This product is for in vitro research use only