

Dynein light chain 1, cytoplasmic

Catalog No: #AP82804



Package Size: #AP82804-1 50ug #AP82804-2 100ug #AP82804-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dynein light chain 1, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Dynein light chain LC8-type 1, DNCL1
Accession No.	Uniprot ID:P61285
Uniprot	P61285
GeneID	101902172;404151;784058;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DYNLL1

Note: This product is for in vitro research use only