

Coagulation factor XIII B chain

Catalog No: #AP82805



Package Size: #AP82805-1 50ug #AP82805-2 100ug #AP82805-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Coagulation factor XIII B chain |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Protein-glutamine gamma-glutamyltransferase B chain,Transglutaminase B chain,Cf13b |
| Accession No. | Uniprot ID:Q07968 |
| Uniprot | Q07968 |
| GeneID | 14060; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:F13b

Note: This product is for in vitro research use only