

Inactive dipeptidyl peptidase 10

Catalog No: #AP82806



Package Size: #AP82806-1 50ug #AP82806-2 100ug #AP82806-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Inactive dipeptidyl peptidase 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Dipeptidyl peptidase X,DPP X,Kv4 potassium channel auxiliary subunit
Accession No.	Uniprot ID:Q6Q629
Uniprot	Q6Q629
GeneID	363972;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dpp10

Note: This product is for in vitro research use only