

D(2) dopamine receptor

Catalog No: #AP82808



Package Size: #AP82808-1 50ug #AP82808-2 100ug #AP82808-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	D(2) dopamine receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Dopamine D2 receptor
Accession No.	Uniprot ID:P20288
Uniprot	P20288
GeneID	281126;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DRD2

Note: This product is for in vitro research use only