

Protein ripply3

Catalog No: #AP82809



Package Size: #AP82809-1 50ug #AP82809-2 100ug #AP82809-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein ripply3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Down syndrome critical region protein 6,Dscr6
Accession No.	Uniprot ID:Q924S9
Uniprot	Q924S9
GeneID	170765;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ripply3

Note: This product is for in vitro research use only