

Ephrin type-B receptor 4

Catalog No: #AP82810



Package Size: #AP82810-1 50ug #AP82810-2 100ug #AP82810-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Ephrin type-B receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Developmental kinase 2,mDK-2,Hepatoma transmembrane kinase,Tyrosine kinase MYK-1,Htk,Mdk2,Myk1
Accession No.	Uniprot ID:P54761
Uniprot	P54761
GeneID	13846;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ephb4

Note: This product is for in vitro research use only