

Dynein light chain 1, cytoplasmic

Catalog No: #AP82813



Package Size: #AP82813-1 50ug #AP82813-2 100ug #AP82813-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dynein light chain 1, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	8 kDa dynein light chain,DLC8,Dynein light chain LC8-type 1,Protein inhibitor of neuronal nitric oxide synthase,PIN,Dlc1,Dncl1,Dncl1
Accession No.	Uniprot ID:P63168
Uniprot	P63168
GeneID	56455;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dynll1

Note: This product is for in vitro research use only