

Dynamin-2

Catalog No: #AP82814



Package Size: #AP82814-1 50ug #AP82814-2 100ug #AP82814-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dynamin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dynamin UDNM,Dyn2
Accession No.	Uniprot ID:P39054
Uniprot	P39054
GeneID	13430;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dnm2

Note: This product is for in vitro research use only