

# EGF-like repeat and discoidin I-like domain-containing protein 3



Catalog No: #AP82815

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP82815-1 50ug #AP82815-2 100ug #AP82815-3 1mg

## Description

Product Name	EGF-like repeat and discoidin I-like domain-containing protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Developmentally-regulated endothelial cell locus 1 protein,Integrin-binding protein DEL1,Del1
Accession No.	Uniprot ID:O35474
Uniprot	O35474
GeneID	13612;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Edil3

Note: This product is for in vitro research use only