

## Eukaryotic translation initiation factor 1

Catalog No: #AP82816



Package Size: #AP82816-1 50ug #AP82816-2 100ug #AP82816-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Eukaryotic translation initiation factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	eIF1, Protein translation factor SUI1 homolog, Sui1, Sui1-rs1
Accession No.	Uniprot ID:P48024
Uniprot	P48024
GeneID	20918;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Eif1

Note: This product is for in vitro research use only