

## Fibronectin type III domain-containing protein 3B

Catalog No: #AP82818



Package Size: #AP82818-1 50ug #AP82818-2 100ug #AP82818-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fibronectin type III domain-containing protein 3B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Factor for adipocyte differentiation 104,HCV NS5A-binding protein 37,Fad104,Kiaa4164,Ns5abp37
Accession No.	Uniprot ID:Q6NWW9
Uniprot	Q6NWW9
GeneID	72007;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Fndc3b

Note: This product is for in vitro research use only