

Ectonucleoside triphosphate diphosphohydrolase 5

Catalog No: #AP82819



Package Size: #AP82819-1 50ug #AP82819-2 100ug #AP82819-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ectonucleoside triphosphate diphosphohydrolase 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	NTPDase 5, Guanosine-diphosphatase ENTPD5, GDPase ENTPD5, Uridine-diphosphatase ENTPD5, UDPase ENTPD5
Accession No.	Uniprot ID:Q6P6S9
Uniprot	Q6P6S9
GeneID	314312;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Entpd5

Note: This product is for in vitro research use only