

Early endosome antigen 1

Catalog No: #AP82824



Package Size: #AP82824-1 50ug #AP82824-2 100ug #AP82824-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Early endosome antigen 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Endosome-associated protein p162,Zinc finger FYVE domain-containing protein 2,ZFYVE2
Accession No.	Uniprot ID:Q15075
Uniprot	Q15075
GeneID	8411;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EEA1

Note: This product is for in vitro research use only