

Eukaryotic elongation factor 2 kinase

Catalog No: #AP82826



Package Size: #AP82826-1 50ug #AP82826-2 100ug #AP82826-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Eukaryotic elongation factor 2 kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	eEF-2 kinase, Calcium/calmodulin-dependent eukaryotic elongation factor 2 kinase
Accession No.	Uniprot ID:O08796
Uniprot	O08796
GeneID	13631;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Eef2k

Note: This product is for in vitro research use only