

## Endoplasmin

Catalog No: #AP82828



Package Size: #AP82828-1 50ug #AP82828-2 100ug #AP82828-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Endoplasmin  |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Pig  |
| Other Names        | 94 kDa glucose-regulated protein,GRP-94,98 kDa protein kinase,PPK 98,Heat shock protein 90 kDa beta member 1,gp96 homolog,GRP94,TRA1 |
| Accession No.      | Uniprot ID:Q29092  |
| Uniprot            | Q29092   |
| GeneID             | 397191;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)   |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HSP90B1

Note: This product is for in vitro research use only