

## Ephrin type-A receptor 4

Catalog No: #AP82829



Package Size: #AP82829-1 50ug #AP82829-2 100ug #AP82829-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ephrin type-A receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EPH-like kinase 8,EK8,Tyrosine-protein kinase TYRO1,Tyrosine-protein kinase receptor SEK,HEK8,SEK,TYRO1
Accession No.	Uniprot ID:P54764
Uniprot	P54764
GeneID	2043;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:EPHA4

Note: This product is for in vitro research use only