

Fibroblast growth factor 7

Catalog No: #AP82833



Package Size: #AP82833-1 50ug #AP82833-2 100ug #AP82833-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FGF-7,Heparin-binding growth factor 7,HBGF-7,Keratinocyte growth factor,KGF,Fgf-7,Kgf
Accession No.	Uniprot ID:P36363
Uniprot	P36363
GeneID	14178;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fgf7

Note: This product is for in vitro research use only