

Enteropeptidase

Catalog No: #AP82835



Package Size: #AP82835-1 50ug #AP82835-2 100ug #AP82835-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Enteropeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Enterokinase,Serine protease 7,Transmembrane protease serine 15,ENTK,PRSS7
Accession No.	Uniprot ID:P98074
Uniprot	P98074
GeneID	397152;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TMPRSS15

Note: This product is for in vitro research use only