

Neutrophil elastase

Catalog No: #AP82840



Package Size: #AP82840-1 50ug #AP82840-2 100ug #AP82840-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neutrophil elastase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Elastase-2,Leukocyte elastase,Ela2
Accession No.	Uniprot ID:Q3UP87
Uniprot	Q3UP87
GeneID	50701;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Elane

Note: This product is for in vitro research use only