

## Eukaryotic translation initiation factor 2-alpha kinase 3

Catalog No: #AP82841



Package Size: #AP82841-1 50ug #AP82841-2 100ug #AP82841-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Eukaryotic translation initiation factor 2-alpha kinase 3                    |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Ms   |
| Other Names        | PRKR-like endoplasmic reticulum kinase,Pancreatic eIF2-alpha kinase,Pek,Perk |
| Accession No.      | Uniprot ID:Q9Z2B5  |
| Uniprot            | Q9Z2B5   |
| GeneID             | 13666;   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                       |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Eif2ak3

Note: This product is for in vitro research use only