

Ephrin-A1

Catalog No: #AP82842



Package Size: #AP82842-1 50ug #AP82842-2 100ug #AP82842-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ephrin-A1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	EPH-related receptor tyrosine kinase ligand 1, LERK-1, Immediate early response protein B61, Epgl1, Ep11, Lerk1
Accession No.	Uniprot ID: P52793
Uniprot	P52793
GeneID	13636;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Efn1

Note: This product is for in vitro research use only