

Fatty acid-binding protein, heart

Catalog No: #AP82844



Package Size: #AP82844-1 50ug #AP82844-2 100ug #AP82844-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fatty acid-binding protein, heart
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Fatty acid-binding protein 3,Heart-type fatty acid-binding protein,H-FABP,Mammary-derived growth inhibitor,MDGI
Accession No.	Uniprot ID:P10790
Uniprot	P10790
GeneID	281758;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FABP3

Note: This product is for in vitro research use only