

Ephrin type-A receptor 3

Catalog No: #AP82846



Package Size: #AP82846-1 50ug #AP82846-2 100ug #AP82846-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ephrin type-A receptor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	EPH-like kinase 4,EK4,Tyrosine-protein kinase receptor ETK1,CEK4,ETK1
Accession No.	Uniprot ID:P29318
Uniprot	P29318
GeneID	396402;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EPHA3

Note: This product is for in vitro research use only