

Sorbitol dehydrogenase

Catalog No: #AP82850



Package Size: #AP82850-1 50ug #AP82850-2 100ug #AP82850-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sorbitol dehydrogenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	L-iditol 2-dehydrogenase
Accession No.	Uniprot ID:Q00796
Uniprot	Q00796
GeneID	6652;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SORD

Note: This product is for in vitro research use only