

DNA ligase 1

Catalog No: #AP82852



Package Size: #AP82852-1 50ug #AP82852-2 100ug #AP82852-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA ligase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DNA ligase I, Polydeoxyribonucleotide synthase [ATP] 1, Lig-1
Accession No.	Uniprot ID:P37913
Uniprot	P37913
GeneID	16881;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Lig1

Note: This product is for in vitro research use only