

Dysbindin

Catalog No: #AP82853



Package Size: #AP82853-1 50ug #AP82853-2 100ug #AP82853-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Dysbindin |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Chicken |
| Other Names | Biogenesis of lysosome-related organelles complex 1 subunit 8,BLOC-1 subunit 8,Dysbindin-1,Dystrobrevin-binding protein 1,Bloc1s8,RCJMB04_10b8 |
| Accession No. | Uniprot ID:Q5ZKM0 |
| Uniprot | Q5ZKM0 |
| GeneID | 420840; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DTNBP1

Note: This product is for in vitro research use only